

Accelerating bioprocess workflows with microfluidic immunoassay-based parallel impurity analysis

Rob Durham, Pirjo Lehtonen and Marie Andersson

Gyros Protein Technologies AB, Uppsala Science Park, Dag Hammarskjölds väg 54, 751 83 Uppsala, Sweden

BACKGROUND

Product titer and impurities levels are critical quality attributes (CQA) in biopharmaceutical manufacturing. Minimizing the level of bioprocess impurities in the downstream purification process is a critical regulatory and safety requirement since they can induce adverse effects in patients. As the number of impurities requiring analysis continues to expand, executing accurate and precise bioassays of different types for large sample sets can cause workflow bottlenecks.

The Gyrolab® xPand system enables methods to be mixed and matched to meet specific needs, either by performing up to five CD runs with one assay and many samples, or by setting up a Gyroplex® panel, when up to 5 independent immunoassays are run sequentially in an automated walk-away format, reducing hands-on time and increasing productivity in bioprocess analytics. Since each assay is run independently with no cross-talk, re-runs of a single analyte is easy without having to sacrifice reagents and samples to repeat the assay for the other analytes.

RESULTS

In the present study, a MabSelect SuRe™ purification eluate from a Chinese Hamster Ovary (CHO) mAb culture was assessed for IgG titer and four different process related impurity levels in a Gyroplex panel generating results in 5 hours. (Data show dilution linearity and spike recovery, which are important assay validation parameters in all stages of the biopharmaceutical workflow (Fig. 2 and Table 1).) The Gyrolab software automatically calculates spike recovery from the dilution series of spiked and unspiked unknowns for each assay. This enables the accuracy and dilution linearity of sample analysis to be verified. The system helps the user to set criteria for dilutions and spikes according to current guidelines (e.g. USP 38-NF 33 1132 Residual HCP; see Fig. 1).





Figure 2. Dilution linearity plots (A), spike recovery and dilution linearity data (Conc vs Max (%) column) (B) and Standard curves (C), for CHO host cell proteins (HCP), Insulin, human Transferrin, Protein A and human IgG.

Table 1. Summary table of impurity levels in eluate after MabSelect SuRe purification

Impurity	Level (ng/mL)	CV (%)	ng impurity/mg IgG product
CHO HCP	10000	8.8	1279
Insulin	2.1	9.8	0.268
Transferrin	0.01	3.8	0.0013
Protein A	52	4.4	6.65

THE GYROLAB SYSTEM

Gyrolab systems perform automated immunoassays within nanoliter-scale microfluidic structures in a Compact Disk (CD) format. Each structure on the CD comprises a 15-nanoliter affinity column pre-packed with streptavidin-coated particles, supporting a variety of assay types including sandwich and indirect antibody assays. While Gyrolab xPlore™ runs single CDs, Gyrolab xPand can run up to five CDs unattended. Consumption of sample and

reagents is dramatically reduced compared with plate-based ELISA. Microfluidic control ensures that all samples on a CD are processed in parallel, giving consistent results. Each microstructure equates to one data point, eliminating cross talk. The control and analysis software enables 21 CFR part 11 compliance, ensuring that assays can be developed and transferred through to GMP and GLP environments.



MATERIALS AND METHODS

Automated three-step sandwich assays were used to measure titer as well as Protein A, CHO HCP, Insulin and transferrin impurities in a MabSelect SuRe™ purified therapeutic antibody sample* from a CHO transfected cell line. The ready-to-use Gyrolab CHO-HCP Kit, Gyrolab Protein A Kit and Gyrolab hulgG – High titer Kit were used. Reagents, CDs and buffers needed for all assays are summarized in Table 2.

Three different types of Gyrolab Bioaffy™ CDs were used, including a Gyrolab Mixing CD 96 in the Protein A kit that includes an automated mixing step to perform acid dissociation prior to measurement of Protein A.

Gyrolab Bioaffy CD




Table 2. Assay designs and assay performances. The hulgG, CHO HCP and Protein A are available as ready-to-use Gyrolab kits and the Insulin and Transferrin assays are available as Gyrolab assay protocols **. The assay performances of each assay are from respective IFU or protocol instruction.

Assay	Capture reagent	Detection reagent	Gyrolab Bioaffy CD	Analytical range (ng/mL)		Precision (CV %)	
				LLOQ	ULOQ	Intra-run	Inter-run
CHO HCP	Affinity purified goat antibody	Affinity purified goat antibody	1000 HC	~2	8 000	< 10	< 20
Insulin	mAb mouse anti insulin/proinsulin #1	mAb mouse anti insulin/proinsulin #2	1000	0.025	15	< 10	< 20
Transferrin (human)	anti-human transferrin mAb	anti-human transferrin pAb	1000	0.1	150	< 10	< 10
Protein A (MabSelect SuRe)	Affinity purified a-Pro A chicken pAb	Affinity purified a-Pro A chicken pAb	Mixing CD 96	~0.1	100	< 10	< 10
IgG Titer	Protein A derivative from Staphylococcus aureus	F(ab) ₂ fragment of anti-human IgG	20 HC	300	1 000 000	< 10	< 10

*Gyros Protein Technologies wishes to thank POLYMUN Scientific Immunobiologische Forschung GmbH for generously supplying the samples for analysis.

** Gyrolab assay protocols for Insulin and human Transferrin can be downloaded from <https://www.gyrosproteintechnologies.com/gyrolab-assays>

Gyrolab, Rexpix and Gyroplex are registered trademarks and Gyros, Gyrolab xPlore, Bioaffy and Gyros logo are trademarks of Gyros Protein Technologies Group. All other trademarks are the property of their respective owners. Products and technologies from Gyros Protein Technologies are covered by one or more patents and/or proprietary intellectual property rights. All infringements are prohibited and will be prosecuted. Please contact Gyros Protein Technologies AB for further details. Products are for research use only. Not for use in diagnostic procedures. © Gyros Protein Technologies AB 2021. D0039913/A